



SURGI EDGETM



SURGICAL BLADES



Complete selection of the most widely used blade size in both - Carton Steel and Stainless Steel. Individually VCI Paper and foil packed.

DISPOSABLE SCALPELS



Pre-assembled Stainless Steel blades mounted on plastic handle with protection guard. Individually pouch packed. Available in different color handle.

SAFETY SCALPELS



Designed to prevent from accidental needle sticks. Blades extends and retracts with just one hand. Ribbed handle for better grip and control with gloved hands. Glides easily and comfortably in to locked position. Metric ruler assists in determining length of incision. Disposable, sterile and ready to use, Available in different color handle.

SKIN GRAFT BLADES



Designed to fit Braithwaite, Cobbett and Waston handles. Individually VCI Paper wrapped and pouch packed.

POSTMORTEM BLADES



Available with numbers 40, 40B, 60, 60B and 62 in Carbon Steel, Individually VCI Paper packed.

BLOOD LANCET



Designed for easy puncturing to minimize patient trauma, Individually packed.

STITCH CUTTERS



Long and medium Stitch Cutter does not require handle. Mini Stitch Cutter fits on plastic handle. Also available pre-assembled Stitch cutters with Plastic handles.

BIOPSY PUNCH



Made from Stainless Steel with sharp cutting edge, Ribbed handle for better grip and control. Size embossed on each punch for easy identification. Disposable, Sterile and ready to use. Available in 1.0 mm, 1.5 mm, 2.0 mm, 3.0 mm, 3.5 mm, 4.0 mm, 5.0 mm, 6.0 mm, 8.0 mm.



SURGI EDGE (India)

305, "VEDANT", Near Ganesh Plaza,
Navrangpura, Ahmedabad - 380 009. INDIA
Phone : +91-79-2646 2104 Fax : +91-79-2642 6847
E-mail : info@surgiedge.com
Website : WWW.surgiedge.com

